

NF-KB p100 Rabbit mAb


ABZ2272

Product Overview


Name	NF-KB p100 Rabbit mAb
Accession(Primary)	Q00653
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 97 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; IF-1:10-1:100

Target information (Q00653)

Gene Name	NFKB2
About Gene Symbol	NFKB2 (Nuclear Factor Kappa B Subunit 2) is a transcription factor involved in NF- κ B signaling and adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic ta



Western blot analysis of NF-KB p100/NF-KB2 in C6, 3T3, HeLa lysates using NF-KB p100 antibody.



Immunocytochemistry analysis of NF-KB p100(green) in HeLa using NF-KB p100 antibody, and DAPI(blue)