

## CDK16 Rabbit mAb


ABZ2273

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | CDK16 Rabbit mAb                           |
| <b>Accession(Primary)</b>   | Q00536                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 56 kDa; Observed MW: 65 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000                           |

### Target information (Q00536)

|                          |                                                                                                                                                                                                                                                                                                   |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | CDK16                                                                                                                                                                                                                                                                                             |
| <b>About Gene Symbol</b> | CDK16 (Cyclin Dependent Kinase 16) is a serine/threonine protein kinase involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of PCTAIRE1 in A549, HL-60, U251 lysates using CDK16 antibody.