

TFIIB Rabbit mAb

ABZ2275


Product Overview


Name	TFIIB Rabbit mAb
Accession(Primary)	Q00403
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q00403)

Gene Name GTF2B

About Gene Symbol





Western blot analysis of TFIIB in HeLa, A549, HL-60, U2OS, C6 lysates using TFIIB antibody.