

HDLBP Rabbit mAb

ABZ2276


Product Overview

Name	HDLBP Rabbit mAb
Accession(Primary)	Q00341
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 141 kDa; Observed MW: 141 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:100-1:200

Target information (Q00341)

Gene Name HDLBP

About Gene Symbol



Western blot analysis of HDLBP in HeLa, A549, HL-60, U251, U87-MG lysates using HDLBP antibody.

