


MAT1A Rabbit mAb

ABZ2277


Product Overview

Name	MAT1A Rabbit mAb
Accession(Primary)	Q00266
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (Q00266)

Gene Name MAT1A**About Gene Symbol**

Western blot analysis of MAT1A in rat liver lysates using MAT1A antibody.



Western blot analysis of MAT1A in HSC, mouse liver lysates using MAT1A antibody.