

mtTFA Rabbit mAb


ABZ2278

Product Overview


Name	mtTFA Rabbit mAb
Accession(Primary)	Q00059
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 29 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; IP-1:10-1:100

Target information (Q00059)


Gene Name	TFAM
About Gene Symbol	



Western blot analysis of mtTFA in K562, HeLa lysates using mtTFA antibody.



Immunocytochemistry analysis of mtTFA/TFAM(green) in HeLa using mtTFA/TFAM antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using mtTFA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.