

Cytochrome C Rabbit mAb


ABZ2279

Product Overview



Name	Cytochrome C Rabbit mAb
Accession(Primary)	P99999
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 12 kDa; Observed MW: 12 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200; ICC/IF-1:50-1:100

Target information (P99999)

Gene Name	CYCS
About Gene Symbol	



Immunocytochemistry analysis of Cytochrome C (green) in K562 using Cytochrome C antibody, and DAPI (blue).



Western blot analysis of Cytochrome C in K562, rat Brain, C6, 3T3, HeLa lysates using Cytochrome C antibody.