

Cytochrome C Rabbit mAb


ABZ2280

Product Overview


Name	Cytochrome C Rabbit mAb
Accession(Primary)	P99999
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 12 kDa; Observed MW: 12 kDa
Recommended Dilution	WB-1:1000-1:10000; IHC-1:200-1:5000; ICC/IF-1:50-1:100; IP-1:20-1:50

Target information (P99999)


Gene Name	CYCS
About Gene Symbol	



Western blot analysis of Cytochrome C in mouse brain lysates using Cytochrome C antibody.



Western blot analysis of Cytochrome C in Jurkat, rat Brain, C6 lysates using Cytochrome C antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Cytochrome C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.