

RBM3 Rabbit mAb


ABZ2281

Product Overview


Name	RBM3 Rabbit mAb
Accession(Primary)	P98179
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100; FC-1:20

Target information (P98179)


Gene Name	RBM3
About Gene Symbol	



Western blot analysis of RBM3 in 3T3, mouse brain lysates using RBM3 antibody.



Western blot analysis of RBM3 in A549, HL-60, U251 lysates using RBM3 antibody.



Immunocytochemistry analysis of RBM3 (green) in
hela using RBM3 antibody, and DAPI(blue)