

Caspase 7 Rabbit mAb


ABZ2283

Product Overview

Name	Caspase 7 Rabbit mAb
Accession(Primary)	P97864
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; IP-1:50-1:100

Target information (P97864)

Gene Name	Casp7
About Gene Symbol	



Western blot analysis of Caspase7 in C6, 3T3
lysates using Caspase 7 antibody.