

SRSF3 Rabbit mAb


ABZ2285

Product Overview


Name	SRSF3 Rabbit mAb
Accession(Primary)	P84103
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 19 kDa; Observed MW: 19 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:100-1:200

Target information (P84103)


Gene Name	SRSF3
About Gene Symbol	



Western blot analysis of SFRS3 in mouse liver, mouse spleen lysates using SRSF3 antibody.



Western blot analysis of SFRS3 in Jurkat, C6, HeLa lysates using SFRS3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using SFRS3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.