

Smad3 Rabbit mAb


ABZ2287

Product Overview


Name	Smad3 Rabbit mAb
Accession(Primary)	P84022
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Hamster
Calculated MW	Calculated MW: 48 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:2000-1:4000; IHC-1:50-1:100; IP-1:20

Target information (P84022)


Gene Name	SMAD3
About Gene Symbol	SMAD3 (SMAD Family Member 3) is a transcription factor involved in TGF- β signaling. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Smad3 in 3T3, C6 lysates using Smad3 antibody.



Western blot analysis of Smad3 in Jurkat lysates using Smad3 antibody.



Western blot analysis of Smad3 in K562, CHO-K1
lysates using Smad3 antibody