

ABAT Rabbit mAb


ABZ2288

Product Overview

Name	ABAT Rabbit mAb
Accession(Primary)	P80404
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 56 kDa; Observed MW: 56 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; IP-1:10-1:100

Target information (P80404)

Gene Name	ABAT
About Gene Symbol	



Western blot analysis of ABAT/GABA-T in rat Brain lysates using ABAT antibody.