

ABAT Rabbit mAb


ABZ2289

Product Overview


Name	ABAT Rabbit mAb
Accession(Primary)	P80404
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 56 kDa; Observed MW: 56 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1; 50-1:100; FC-1:50-1:100; IP-1:20-1:50

Target information (P80404)


Gene Name	ABAT
About Gene Symbol	



Western blot analysis of ABAT/GABA-T in MCF-7, mouse liver lysates using ABAT antibody.



Western blot analysis of GABA-T in rat Brain lysates using GABA-T antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using ABAT/GABAT antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.