


MS2 Rabbit mAb


ABZ2294

Product Overview

Name	MS2 Rabbit mAb
Accession(Primary)	P78325
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 89 kDa; Observed MW: 89 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:50

Target information (P78325)

Gene Name ADAM8**About Gene Symbol**



Western blot analysis of MS2 in K562 lysates
using MS2 antibody.