

SIRP alpha Rabbit mAb


ABZ2295

Product Overview


Name	SIRP alpha Rabbit mAb
Accession(Primary)	P78324
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 85 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:200-1:500; FC-1:20-1:50

Target information (P78324)

Gene Name	SIRPA
About Gene Symbol	



Western blot analysis of SIRP alpha in mouse brain lysates using SIRP alpha antibody.



Western blot analysis of SIRP alpha in rat Brain lysates using SIRP alpha antibody.