

SIRP alpha Rabbit mAb


ABZ2296

Product Overview


Name	SIRP alpha Rabbit mAb
Accession(Primary)	P78324
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 55 kDa; Observed MW: 85 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200; FC-1:200-1:500

Target information (P78324)

Gene Name	SIRPA
About Gene Symbol	



Western blot analysis of SIRP alpha in A549, MCF-7, 293T lysates using SIRP alpha antibody.



Immunocytochemistry analysis of SIRP alpha (green) in HeLa using SIRP alpha antibody, and DAPI (blue)