

Caspase 11 Rabbit mAb


ABZ2297

Product Overview

Name	Caspase 11 Rabbit mAb
Accession(Primary)	P70343
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43,38,25 kDa
Recommended Dilution	WB-1:2000-1:10000; IP-1:20-1:50

Target information (P70343)

Gene Name	Casp4
About Gene Symbol	



Western blot analysis of Caspase11 in C6, 3T3, Hela lysates using Caspase 11 antibody.