

ACTA1 Rabbit mAb

ABZ2300

Product Overview

Name	ACTA1 Rabbit mAb
Accession(Primary)	P68133
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:50-1:100

Target information (P68133)

Gene Name ACTA1

About Gene Symbol



