

ACTA1 Rabbit mAb


ABZ2301

Product Overview


Name	ACTA1 Rabbit mAb
Accession(Primary)	P68133
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (P68133)


Gene Name	ACTA1
About Gene Symbol	




Western blot analysis of Alpha Skeletal Muscle in rat muscle lysates using alpha Skeletal Muscle Actin antibody.



Western blot analysis of muscle Actin in Hela lysates using muscle Actin antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using muscle Actin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of muscle Actin(green) in Hela using muscle Actin antibody, and DAPI(blue).

