

ACTA1 Rabbit mAb


ABZ2302

Product Overview


Name	ACTA1 Rabbit mAb
Accession(Primary)	P68133
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:20

Target information (P68133)


Gene Name	ACTA1
About Gene Symbol	



Immunocytochemistry analysis of alpha Skeletal Muscle Actin (green) in HeLa using alpha Skeletal Muscle Actin antibody, and DAPI (blue).



Western blot analysis of muscle Actin in HeLa, CHO-K1, C6, rat Brain lysates using muscle Actin antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using muscle Actin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.