

YB1 Rabbit mAb


ABZ2307

Product Overview


Name	YB1 Rabbit mAb
Accession(Primary)	P67809
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 49 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:200; ICC/IF-1:50-1:100; FC-1:50-1:100; IP-1:20-1:50

Target information (P67809)


Gene Name	YBX1
About Gene Symbol	




Immunocytochemistry analysis of YB1 (green) in Hela using YB1 antibody, and DAPI (blue).



Western blot analysis of YB1 in mouse heart lysates using YB1 antibody.



Western blot analysis of YB1 in Jurkat, Hela lysates using YB1 antibody



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using YB1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.