

YB1 Rabbit mAb


ABZ2308

Product Overview


Name	YB1 Rabbit mAb
Accession(Primary)	P67809
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 36 kDa; Observed MW: 49 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (P67809)


Gene Name	YBX1
About Gene Symbol	




Western blot analysis of YB1 in 3T3 lysates using YB1 antibody.



Western blot analysis of YB1 in HeLa, C2C12 lysates using YB1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using YB1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of YB1 (green) in hela using YB1 antibody, and DAPI (blue).