

RACK1 Rabbit mAb


ABZ2311

Product Overview


Name	RACK1 Rabbit mAb
Accession(Primary)	P63244
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000

Target information (P63244)

Gene Name	RACK1
About Gene Symbol	RACK1 (Receptor For Activated C Kinase 1) is a membrane receptor involved in translational regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of RACK1 in C6, 3T3, HeLa lysates using RACK1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using RACK1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.