

eIF5A Rabbit mAb


ABZ2312

Product Overview


Name	eIF5A Rabbit mAb
Accession(Primary)	P63241
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:500; ICC/IF-1:200-1:500; FC-1:50-1:100

Target information (P63241)


Gene Name	EIF5A
About Gene Symbol	



Western blot analysis of eIF5A in Jurkat lysates using eIF5A antibody.



Western blot analysis of eIF5A in Jurkat, THP-1, C6, 3T3, HeLa lysates using eIF5A antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using eIF5A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.