

Skp1 Rabbit mAb


ABZ2313

Product Overview


Name	Skp1 Rabbit mAb
Accession(Primary)	P63208
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 19 kDa; Observed MW: 19 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:200-1:500; FC-1:100-1:1000

Target information (P63208)

Gene Name	SKP1
About Gene Symbol	



Western blot analysis of Skp1 in K562, rat Brain, C6, 3T3 lysates using SKP1 antibody.



Immunocytochemistry analysis of Skp1 (green) in MCF-7 using Skp1 antibody, and DAPI (blue)