

DYNLL1/PIN Rabbit mAb


ABZ2314

Product Overview


Name	DYNLL1/PIN Rabbit mAb
Accession(Primary)	P63167
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 10 kDa; Observed MW: 10 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200; FC-1:200-1:2000; IP-1:20-1:50

Target information (P63167)

Gene Name	DYNLL1
About Gene Symbol	DYNLL1 (Dynein Light Chain LC8-Type 1) is a biosynthetic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of DYNLL1 in 3T3, HeLa lysates using DYNLL1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using DYNLL1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.