

Ubiquitin 60S Ribosomal Protein L40 Rabbit mAb


ABZ2318

Product Overview


Name	Ubiquitin 60S Ribosomal Protein L40 Rabbit mAb
Accession(Primary)	P62987
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 15 kDa; Observed MW: 10 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:200-1:500

Target information (P62987)

Gene Name	UBA52
About Gene Symbol	



Immunocytochemistry analysis of Ubiquitin 60S Ribosomal Protein L40 (green) in HeLa using Ubiquitin 60S Ribosomal Protein L40 antibody, and DAPI (blue).



Western blot analysis of UBA52 in HeLa, A549, HL-60 lysates using UBA52 antibody.