

GNB2 Rabbit mAb


ABZ2319

Product Overview


Name	GNB2 Rabbit mAb
Accession(Primary)	P62879
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-1:50-1:100

Target information (P62879)

Gene Name	GNB2
About Gene Symbol	GNB2 (Transducin Beta Chain 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellu



Western blot analysis of GNB2 in rat Brain, C6, 3T3, HeLa lysates using GNB2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using GNB2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.