

Ran Rabbit mAb


ABZ2321

Product Overview


Name	Ran Rabbit mAb
Accession(Primary)	P62826
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:100-1:200; FC-1:50-1:100; IP-1:20-1:50

Target information (P62826)


Gene Name	RAN
About Gene Symbol	



Immunocytochemistry analysis of Ran (green) in U87-MG using Ran antibody, and DAPI (blue).



Western blot analysis of Ran in Raw264.7 lysates using Ran antibody.



Western blot analysis of Ran in Jurkat, HeLa,
COS7, Raji lysates using Ran antibody