

Histone H4 Rabbit mAb


ABZ2322

Product Overview


Name	Histone H4 Rabbit mAb
Accession(Primary)	P62805
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 11 kDa; Observed MW: 11 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:50-1:100

Target information (P62805)


Gene Name	H4C1
About Gene Symbol	



Immunocytochemistry analysis of Histone H4 (green) in HeLa using Histone H4 antibody, and DAPI (blue).



Western blot analysis of Histone H4 in HeLa lysates using Histone H4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using Histone H4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.