


SNRPD3 Rabbit mAb

ABZ2325

Product Overview

Name	SNRPD3 Rabbit mAb
Accession(Primary)	P62318
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P62318)

Gene Name SNRPD3**About Gene Symbol**

Western blot analysis of SNRPD3 in C6, HeLa lysates using SNRPD3 antibody.

