

RPS8 Rabbit mAb


ABZ2326

Product Overview


Name	RPS8 Rabbit mAb
Accession(Primary)	P62241
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:200-1:500; FC-1:100-1:200; IP-1:20-1:50

Target information (P62241)


Gene Name	RPS8
About Gene Symbol	



Western blot analysis of RPS8 in Raw264.7, PC-12 lysates using RPS8 antibody.



Western blot analysis of RPS8 in HeLa, CHO-K1, C6, Jurkat lysates using RPS8 antibody.



Immunocytochemistry analysis of RPS8 (green) in
hela using RPS8 antibody, and DAPI (blue)