

PSMC5 Rabbit mAb


ABZ2327

Product Overview


Name	PSMC5 Rabbit mAb
Accession(Primary)	P62195
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:100-1:500; FC-1:100-1:500; IP-1:10-1:100

Target information (P62195)


Gene Name	PSMC5
About Gene Symbol	



Western blot analysis of PSMC5 in mouse heart lysates using PSMC5 antibody.



Western blot analysis of PSMC5 in C6, HeLa lysates using PSMC5 antibody.



Immunocytochemistry analysis of PSMC5(green)
in HeLa using PSMC5 antibody, and DAPI(blue)