

PPP1CB Rabbit mAb


ABZ2328

Product Overview


Name	PPP1CB Rabbit mAb
Accession(Primary)	P62140
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 37 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:50-1:100

Target information (P62140)


Gene Name	PPP1CB
About Gene Symbol	



Western blot analysis of PPP1CB in K562, rat Brain, C6, 3T3, HeLa lysates using PP1C beta antibody.



Immunocytochemistry analysis of PPP1CB(green) in HeLa using PPP1CB antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using PPP1CB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.