

PP1C alpha/beta Rabbit mAb


ABZ2329

Product Overview


Name	PP1C alpha/beta Rabbit mAb
Accession(Primary)	P62136
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:200; FC-1:100-1:200; IP-1:20-1:100

Target information (P62136)


Gene Name	PPP1CA
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded mouse colon using PPP1C alpha/beta antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of PPP1A/PPP1CA+PPP1CB in mouse kidney, Raw264.7 lysates using PPP1A/PPP1CA+PPP1CB antibody.



Western blot analysis of PPP1CA+PPP1CB in K562,
rat Brain lysates using PPP1CA+PPP1CB antibody