

hnRNP K Rabbit mAb


ABZ2331

Product Overview


Name	hnRNP K Rabbit mAb
Accession(Primary)	P61978
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 51 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:2000; ICC/IF-1:20-1:50; FC-1:20; IP-1:20

Target information (P61978)


Gene Name	HNRNPK
About Gene Symbol	HNRNPK (Heterogeneous Nuclear Ribonucleoprotein K) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of hnRNP K in 3T3, MEF lysates using hnRNP K antibody.



Western blot analysis of hnRNP K in C6, CHO-K1, HeLa lysates using hnRNP K antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using hnRNP K antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of hnRNP K (green) in MCF-7 using hnRNP K antibody, and DAPI (blue).

