

LMO4 Rabbit mAb

ABZ2332


Product Overview

Name	LMO4 Rabbit mAb
Accession(Primary)	P61968
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 18 kDa; Observed MW: 18 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:100

Target information (P61968)


Gene Name LMO4

About Gene Symbol



Western blot analysis of LMO4 in 293T, Jurkat lysates using LMO4 antibody.

Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of LMO4 in Jurkat lysates
using LMO4 antibody.