

WDR5 Rabbit mAb


ABZ2333

Product Overview

Name	WDR5 Rabbit mAb
Accession(Primary)	P61964
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 37 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; FC-1:10-1:100

Target information (P61964)

Gene Name	WDR5
About Gene Symbol	



Western blot analysis of WDR5 in K562, C6, 3T3, HeLa lysates using WDR5 antibody.