

beta 2 Microglobulin Rabbit mAb


ABZ2336

Product Overview


Name	beta 2 Microglobulin Rabbit mAb
Accession(Primary)	P61769
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; FC-1:20-1:50; IP-1:20-1:50

Target information (P61769)

Gene Name	B2M
About Gene Symbol	B2M (Beta-2-Microglobulin) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunohistochemistry analysis of paraffin-embedded mouse colon using beta 2 Microglobulin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of beta 2 Microglobulin in HeLa lysates using beta 2 Microglobulin antibody.