



Lysozyme Rabbit mAb

ABZ2338


Product Overview

Name	Lysozyme Rabbit mAb
Accession(Primary)	P61626
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200


Target information (P61626)

Gene Name LYZ**About Gene Symbol**

Immunocytochemistry analysis of Lysozyme (green) in HEPG2 using Lysozyme antibody, and DAPI(blue).



Western blot analysis of Lysozyme in HL-60 lysates using Lysozyme antibody.



Western blot analysis of Lysozyme in Jurkat, C6, 3T3 lysates using Lysozyme antibody