

ATP6V0D1/P39 Rabbit mAb


ABZ2340

Product Overview


Name	ATP6V0D1/P39 Rabbit mAb
Accession(Primary)	P61421
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:100-1:200; FC-1:50-1:100; IP-1:20-1:50

Target information (P61421)


Gene Name	ATP6V0D1
About Gene Symbol	



Western blot analysis of ATP6V0D1 in K562, rat Brain, C6, 3T3, HeLa lysates using ATP6V0D1 antibody.



Immunocytochemistry analysis of ATP6V0D1(green) in HeLa using ATP6V0D1 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using ATP6V0D1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.