

Islet 1 Rabbit mAb


ABZ2341

Product Overview

Name	Islet 1 Rabbit mAb
Accession(Primary)	P61371
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 39 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:1000-1:5000; IP-1:10-1:100

Target information (P61371)

Gene Name	ISL1
About Gene Symbol	



Western blot analysis of Islet 1 in K562, rat Brain lysates using Islet 1 antibody.