

Arp3 Rabbit mAb


ABZ2343

Product Overview


Name	Arp3 Rabbit mAb
Accession(Primary)	P61158
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100; FC-1:20-1:50

Target information (P61158)


Gene Name	ACTR3
About Gene Symbol	ACTR3 (Actin Related Protein 3) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunohistochemistry analysis of paraffin-embedded mouse colon using Actin Related Protein 3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Arp3 in K562, rat Brain, C6, 3T3, HeLa lysates using Arp3 antibody.



Immunocytochemistry analysis of Arp3(green) in
Hela cells using Arp3 antibody, and DAPI(blue)