

Ubiquitin Conjugating Enzyme E2 N Rabbit mAb


ABZ2344

Product Overview


Name	Ubiquitin Conjugating Enzyme E2 N Rabbit mAb
Accession(Primary)	P61088
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100

Target information (P61088)

Gene Name	UBE2N
About Gene Symbol	



Immunocytochemistry analysis of Ubiquitin Conjugating Enzyme E2 N (green) in HeLa using Ubiquitin Conjugating Enzyme E2 N antibody, and DAPI (blue).



Western blot analysis of Ube2N/Ubc13 in HL-60, C6 lysates using Ube2N/Ubc13 antibody.