

Ubiquitin Conjugating Enzyme E2 D3 Rabbit mAb

ABZ2345


Product Overview

Name	Ubiquitin Conjugating Enzyme E2 D3 Rabbit mAb
Accession(Primary)	P61077
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (P61077)

Gene Name UBE2D3

About Gene Symbol



Western blot analysis of UBE2D3 in K562, rat Brain, C6, 3T3, HeLa lysates using Ubiquitin Conjugating Enzyme E2 D3 antibody.

