

RAB10 Rabbit mAb


ABZ2346

Product Overview


Name	RAB10 Rabbit mAb
Accession(Primary)	P61026
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:200-1:500; FC-1:200-1:500; IP-1:20-1:50

Target information (P61026)

Gene Name	RAB10
About Gene Symbol	RAB10 (RAB10, Member RAS Oncogene Family) is a small GTPase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of Rab10 (green) in A549 using Rab10 antibody, and DAPI (blue).



Western blot analysis of RAB10 in C6, 3T3, HeLa lysates using RAB10 antibody.