

Rab8A Rabbit mAb


ABZ2347

Product Overview


Name	Rab8A Rabbit mAb
Accession(Primary)	P61006
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500; FC-1:100-1:200

Target information (P61006)


Gene Name	RAB8A
About Gene Symbol	



Western blot analysis of RAB8A in K562, C6, 3T3, HeLa lysates using Rab8A antibody.



Immunocytochemistry analysis of Rab8A(green) in HeLa using Rab8A antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using RAB8A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.