

Proteasome subunit alpha type 6 Rabbit mAb


ABZ2348

Product Overview


Name	Proteasome subunit alpha type 6 Rabbit mAb
Accession(Primary)	P60900
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 27 kDa; Observed MW: 27 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100

Target information (P60900)

Gene Name	PSMA6
About Gene Symbol	PSMA6 (Proteasome 20S Subunit Alpha 6) is a protease enzyme involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and



Western blot analysis of Proteasome subunit alpha type 6 in K562, C6, 3T3, HeLa lysates using Proteasome alpha 6 antibody.



Immunocytochemistry analysis of Proteasome subunit alpha type 6 (green) in HeLa using Proteasome subunit alpha type 6 antibody, and DAPI(blue)