

SNAP25 Rabbit mAb


ABZ2349

Product Overview


Name	SNAP25 Rabbit mAb
Accession(Primary)	P60880
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (P60880)


Gene Name	SNAP25
About Gene Symbol	SNAP25 (Synaptosome Associated Protein 25) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of SNAP25 (green) in HeLa using SNAP25 antibody, and DAPI (blue).



Western blot analysis of SNAP25 in mouse brain lysates using SNAP25 antibody.



Western blot analysis of SNAP25 in HL-60, C6
lysates using SNAP25 antibody